

CIRBP Polyclonal Antibody

Cat No: HR1AP11908

For research use only

Overview

Product Name	CIRBP Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	CIRBP A18HNRNP CIRP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Seems to play an essential role in cold-induced suppression of cell proliferation.,induction:By cold stress an in response to DNA damage induced by UV irradiation or UV mimetic agents.,similarity:Contains 1 RRM (RNA recognition motif) domain.,tissue specificity:Ubiquitous.,