

AGFG1 Polyclonal Antibody

Cat No: HR1AP11906

For research use only

Overview

Product Name	AGFG1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	61kDa
GeneID?Human?	AGFG1 HRB RAB RIP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	ArfGAP with FG repeats 1(AGFG1) Homo sapiens The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2008],