

ING5 Polyclonal Antibody

Cat No: HR1AP11903

For research use only

Overview

Product Name	ING5 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	26kDa
GeneID?Human?	ING5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>inhibitor of growth family member 5(ING5) Homo sapiens This gene encodes a tumor suppressor protein that inhibits cell growth and induces apoptosis. This protein contains a PHD-type zinc finger. It interacts with tumor suppressor p53 and p300, a component of the histone acetyl transferase complex, suggesting a role in transcriptional regulation. Alternative splicing and the use of multiple promoters and 3' terminal exons results in multiple transcript variants. [provided by RefSeq, Aug 2016],</p>