

# ING2 Polyclonal Antibody

Cat No: HR1AP11902

For research use only

## Overview

Product Name	ING2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	30kDa
GeneID?Human?	ING2 ING1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	inhibitor of growth family member 2(ING2) Homo sapiens This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2014],