

ATRN Polyclonal Antibody

Cat No: HR1AP11901

For research use only

Overview

Product Name	ATRN Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	157kDa
GeneID?Human?	ATRN KIAA0548 MGCA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>attractin(ATRN) Homo sapiens This gene encodes both membrane-bound and secreted protein isoforms. A membrane-bound isoform exhibits sequence similarity with the mouse mahogany protein, a receptor involved in controlling obesity. A secreted isoform is involved in the initial immune cell clustering during inflammatory responses that may regulate the chemotactic activity of chemokines. [provided by RefSeq, Apr 2016],</p>