

CP4AM Polyclonal Antibody

Cat No: HR1AP11895

For research use only

Overview

Product Name	CP4AM Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	CYP4A22
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cytochrome P450 family 4 subfamily A member 22(CYP4A22) Homo sapiens This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This gene is part of a cluster of cytochrome P450 genes on chromosome 1p33. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015].</p>