

# PO3F2 Polyclonal Antibody

Cat No: HR1AP11889

For research use only

## Overview

Product Name	PO3F2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	POU3F2 BRN2 OCT7 OTF7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	POU class 3 homeobox 2(POU3F2) Homo sapiens This gene encodes a member of the POU-III class of neural transcription factors. The encoded protein is involved in neuronal differentiation and enhances the activation of corticotropin-releasing hormone regulated genes. Overexpression of this protein is associated with an increase in the proliferation of melanoma cells. [provided by RefSeq, Mar 2012],