

PKNX1 Polyclonal Antibody

Cat No: HR1AP11888

For research use only

Overview

Product Name	PKNX1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	47kDa
GeneID?Human?	PKNOX1 PREP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the TALE/MEIS homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,tissue specificity:Ubiquitous. Isoform 2 is expressed in all examined tissues except in bone marrow.,