

SMRD1 Polyclonal Antibody

Cat No: HR1AP11882

For research use only

Overview

Product Name	SMRD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	SMARCD1 BAF60A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily d, member 1(SMARCD1) Homo sapiens The protein encoded by this gene is a member of the SWI/SNF family of proteins, whose members display helicase and ATPase activities and which are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein is part of the large ATP-dependent chromatin remodeling complex SNF/SWI and has sequence similarity to the yeast Swp73 protein. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],</p>