

KLF17 Polyclonal Antibody

Cat No: HR1AP11877

For research use only

Overview

Product Name	KLF17 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	KLF17 ZNF393
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function: Binds G/C-rich sites via its zinc fingers and activates transcription from CACCC-box elements. May be a germ cell-specific transcription factor that plays important roles in spermatid differentiation and oocyte development., similarity: Belongs to the Sp1 C2H2-type zinc-finger protein family., similarity: Contains 3 C2H2-type zinc fingers.,