

# KLF14 Polyclonal Antibody

Cat No: HR1AP11876

For research use only

## Overview

Product Name	KLF14 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	KLF14 BTEB5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Kruppel like factor 14(KLF14) Homo sapiens This intronless gene encodes a member of the Kruppel-like family of transcription factors. The encoded protein functions as a transcriptional co-repressor, and is induced by transforming growth factor-beta (TGF-beta) to repress TGF-beta receptor II gene expression. This gene exhibits imprinted expression from the maternal allele in embryonic and extra-embryonic tissues. [provided by RefSeq, Jul 2013],