

CANT1 Polyclonal Antibody

Cat No: HR1AP11874

For research use only

Overview

Product Name	CANT1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	CANT1 SHAPY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	calcium activated nucleotidase 1(CANT1) Homo sapiens This protein encoded by this gene belongs to the apyrase family. It functions as a calcium-dependent nucleotidase with a preference for UDP. Mutations in this gene are associated with Desbuquois dysplasia with hand anomalies. Alternatively spliced transcript variants have been noted for this gene.[provided by RefSeq, Mar 2010],