

C2D1A Polyclonal Antibody

Cat No: HR1AP11872

For research use only

Overview

Product Name	C2D1A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	CC2D1A AKI1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>coiled-coil and C2 domain containing 1A(CC2D1A) Homo sapiens This gene encodes a transcriptional repressor that binds to a conserved 14-bp 5'&#x2013;repressor element and regulates expression of the 5-hydroxytryptamine (serotonin) receptor 1A gene in neuronal cells. The DNA binding and transcriptional repressor activities of the protein are inhibited by calcium. A mutation in this gene results in nonsyndromic mental retardation-3.[provided by RefSeq, Oct 2009],</p>