

BCL7C Polyclonal Antibody

Cat No: HR1AP11871

For research use only

Overview

Product Name	BCL7C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	BCL7C
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	BCL tumor suppressor 7C(BCL7C) Homo sapiens This gene is identified by the similarity of its product to the N-terminal region of BCL7A protein. The BCL7A protein is encoded by the gene known to be directly involved in a three-way gene translocation in a Burkitt lymphoma cell line. The function of this gene has not yet been determined. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2013],