

REV1 Polyclonal Antibody

Cat No: HR1AP11865

For research use only

Overview

Product Name	REV1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	137kDa
GeneID?Human?	REV1 REV1L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>REV1, DNA directed polymerase(REV1) Homo sapiens This gene encodes a protein with similarity to the <i>S. cerevisiae</i> mutagenesis protein Rev1. The Rev1 proteins contain a BRCT domain, which is important in protein-protein interactions. A suggested role for the human Rev1-like protein is as a scaffold that recruits DNA polymerases involved in translesion synthesis (TLS) of damaged DNA. [provided by RefSeq, Mar 2016],</p>