

## CR3L3 Polyclonal Antibody

Cat No: HR1AP11866

For research use only

### Overview

Product Name	CR3L3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	CREB3L3 CREBH HYST1481
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>cAMP responsive element binding protein 3 like 3(CREB3L3) Homo sapiens This gene encodes a member of the basic-leucine zipper family and the AMP-dependent transcription factor family. The encoded protein is localized to the endoplasmic reticulum and acts as a transcription factor activated by cyclic AMP stimulation. The encoded protein binds the cyclic AMP response element (CRE) and the box-B element and has been linked to acute inflammatory response, hepatocellular carcinoma, triglyceride metabolism, and hepcidin expression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012].</p>