

MARE2 Polyclonal Antibody

Cat No: HR1AP11862

For research use only

Overview

Product Name	MARE2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GeneID?Human?	MAPRE2 RP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>microtubule associated protein RP/EB family member 2(MAPRE2) Homo sapiens The protein encoded by this gene shares significant homology to the adenomatous polyposis coli (APC) protein-binding EB1 gene family. This protein is a microtubule-associated protein that is necessary for spindle symmetry during mitosis. It is thought to play a role in the tumorigenesis of colorectal cancers and the proliferative control of normal cells. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jan 2012],</p>