

MCM9 Polyclonal Antibody

Cat No: HR1AP11860

For research use only

Overview

Product Name	MCM9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	125kDa
GeneID?Human?	MCM9 C6orf61 MCMDC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>minichromosome maintenance 9 homologous recombination repair factor(MCM9) Homo sapiens The protein encoded by this gene is a member of the mini-chromosome maintenance (MCM) protein family that are essential for the initiation of eukaryotic genome replication. Binding of this protein to chromatin has been shown to be a pre-requisite for recruiting the MCM2-7 helicase to DNA replication origins. This protein also binds, and is a positive regulator of, the chromatin licensing and DNA replication factor 1, CDT1. [provided by RefSeq, Nov 2010],</p>