

KIF23 Polyclonal Antibody

Cat No: HR1AP11858

For research use only

Overview

Product Name	KIF23 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	105kDa
GeneID?Human?	KIF23 KNSL5 MKLP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>kinesin family member 23(KIF23) Homo sapiens The protein encoded by this gene is a member of kinesin-like protein family. This family includes microtubule-dependent molecular motors that transport organelles within cells and move chromosomes during cell division. This protein has been shown to cross-bridge antiparallel microtubules and drive microtubule movement in vitro. Alternate splicing of this gene results in multiple transcript variants. [provided by RefSeq, Jul 2013],</p>