

CNNM3 Polyclonal Antibody

Cat No: HR1AP11857

For research use only

Overview

Product Name	CNNM3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	77kDa
GeneID?Human?	CNNM3 ACDP3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Probable metal transporter.,miscellaneous:Shares weak sequence similarity with the cyclin family, explaining its name. However it has no cyclin-like function in vivo.,similarity:Belongs to the ACDP family.,similarity:Contains 2 CBS domains.,tissue specificity:Widely expressed. Expressed at higher level in heart and spleen.,