

CNNM1 Polyclonal Antibody

Cat No: HR1AP11856

For research use only

Overview

Product Name	CNNM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	104kDa
GeneID?Human?	CNNM1 ACDP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	cyclin and CBS domain divalent metal cation transport mediator 1(CNNM1) Homo sapiens This gene encodes a member of the ancient conserved domain protein family. The encoded protein may bind copper. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2016],