

NSA2 Polyclonal Antibody

Cat No: HR1AP11851

For research use only

Overview

Product Name	NSA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	28kDa
GeneID?Human?	NSA2 TINP1 HUSSY-29
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>NSA2, ribosome biogenesis homolog(NSA2) Homo sapiens This gene encodes a nucleolar protein involved in cell cycle regulation and proliferation. This gene was identified based on sequence similarity to a highly conserved <i>Saccharomyces cerevisiae</i> gene encoding a pre-ribosomal protein, which is involved in large ribosomal subunit biogenesis. The encoded protein is found at elevated levels in diabetic nephropathy. Alternative splicing results in multiple transcript variants. Several related pseudogenes have been identified. [provided by RefSeq, Nov 2012],</p>