

## XCL2 Polyclonal Antibody

Cat No: HR1AP11852

For research use only

### Overview

Product Name	XCL2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	12kDa
GeneID?Human?	XCL2 SCYC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Chemotactic activity for lymphocytes but not for monocytes or neutrophils.,online information:XCL2 entry,similarity:Belongs to the intercrine gamma family.,