

## TBA8 Polyclonal Antibody

Cat No: HR1AP11849

For research use only

### Overview

Product Name	TBA8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	TUBA8 TUBAL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>TUBA8 encodes a member of the alpha tubulin protein family. Alpha tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. Mutations in this gene are associated with polymicrogyria and optic nerve hypoplasia. Alternate splicing results in multiple transcript variants. TUBA8 (Tubulin Alpha 8) is a Protein Coding gene. Diseases associated with TUBA8 include Cortical Dysplasia, Complex, With Other Brain Malformations 8 and Optic Nerve Hypoplasia. Among its related pathways are Development Slit-Robo signaling and Cytoskeleton remodeling Neurofilaments.</p>