

## TBD Polyclonal Antibody

Cat No: HR1AP11850

For research use only

### Overview

Product Name	TBD Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	TUBD1 TUBD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	alternative products:Experimental confirmation may be lacking for some isoforms,function:In the elongating spermatid it is associated with the manchette, a specialized microtubule system present during reshaping of the sperm head.,similarity:Belongs to the tubulin family.,subcellular location:Associated with centrioles. Both cytoplasmic and nuclear.,