



SP130 Polyclonal Antibody

Cat No: HR1AP11848

For research use only

Overview

Product Name	SP130 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	115kDa
GeneID?Human?	SAP130
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Sin3A associated protein 130(SAP130) Homo sapiens SAP130 is a subunit of the histone deacetylase (see HDAC1; MIM 601241)-dependent SIN3A (MIM 607776) corepressor complex (Fleischer et al., 2003 [PubMed 12724404]).[supplied by OMIM, Mar 2008],