

# MCF2 Polyclonal Antibody

Cat No: HR1AP11841

For research use only

## Overview

Product Name	MCF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GeneID?Human?	MCF2 DBL
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	MCF.2 cell line derived transforming sequence(MCF2) Homo sapiens The oncogenic protein encoded by this gene is a guanine nucleotide exchange factor (GEF) that exerts control over some members of the Rho family of small GTPases. Several transcript variants encoding different isoforms have been found for this gene. These isoforms exhibit different expression patterns and varying levels of GEF activity.[provided by RefSeq, Jan 2010],