

KRT85 Polyclonal Antibody

Cat No: HR1AP11839

For research use only

Overview

Product Name	KRT85 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	KRT85 KRTHB5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>keratin 85(KRT85) Homo sapiens The protein encoded by this gene is a member of the keratin gene family. As a type II hair keratin, it is a basic protein which heterodimerizes with type I keratins to form hair and nails. The type II hair keratins are clustered in a region of chromosome 12q13 and are grouped into two distinct subfamilies based on structure similarity. One subfamily, consisting of KRTHB1, KRTHB3, and KRTHB6, is highly related. The other less-related subfamily includes KRTHB2, KRTHB4, and KRTHB5. [provided by RefSeq, Jul 2008].</p>