

KRT86 Polyclonal Antibody

Cat No: HR1AP11840

For research use only

Overview

Product Name	KRT86 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	53kDa
GeneID?Human?	KRT86 KRTHB6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 86(KRT86) Homo sapiens This gene encodes a type II keratin protein, which heterodimerizes with type I keratins to form hair and nails. This gene is present in a cluster of related genes and pseudogenes on chromosome 12. Mutations in this gene have been observed in patients with the hair disease monilethrix. [provided by RefSeq, Feb 2016],