

K2C6C Polyclonal Antibody

Cat No: HR1AP11836

For research use only

Overview

Product Name	K2C6C Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	62kDa
GeneID?Human?	KRT6C KRT6E
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 6C(KRT6C) Homo sapiens Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. The type II keratins are clustered in a region of chromosome 12q13. [provided by RefSeq, Jul 2009],