

KRT36 Polyclonal Antibody

Cat No: HR1AP11833

For research use only

Overview

Product Name	KRT36 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	51kDa
GeneID?Human?	KRT36 HHA6 HKA6 KRTHA6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	keratin 36(KRT36) Homo sapiens The protein encoded by this gene is a member of the keratin gene family. This type I hair keratin is an acidic protein which heterodimerizes with type II keratins to form hair and nails. The type I hair keratins are clustered in a region of chromosome 17q12-q21 and have the same direction of transcription. [provided by RefSeq, Jul 2008],