

# KDM5D Polyclonal Antibody

Cat No: HR1AP11831

For research use only

## Overview

Product Name	KDM5D Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	169kDa
GeneID?Human?	KDM5D HY HYA JARID1D KIAA0234 SMCY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lysine demethylase 5D(KDM5D) Homo sapiens This gene encodes a protein containing zinc finger domains. A short peptide derived from this protein is a minor histocompatibility antigen which can lead to graft rejection of male donor cells in a female recipient. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009],