

# ITA10 Polyclonal Antibody

Cat No: HR1AP11826

For research use only

## Overview

Product Name	ITA10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	128kDa
GeneID?Human?	ITGA10 UNQ468/PRO827
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>integrin subunit alpha 10(ITGA10) Homo sapiens Integrins are integral transmembrane glycoproteins composed of noncovalently linked alpha and beta chains. They participate in cell adhesion as well as cell-surface mediated signalling. This gene encodes an integrin alpha chain and is expressed at high levels in chondrocytes, where it is transcriptionally regulated by AP-2epsilon and Ets-1. The protein encoded by this gene binds to collagen. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2014],</p>