

# INSL6 Polyclonal Antibody

Cat No: HR1AP11825

For research use only

## Overview

Product Name	INSL6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	23kDa
GeneID?Human?	INSL6 RIF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>insulin like 6(INSL6) Homo sapiens The protein encoded by this gene contains a classical signature of the insulin superfamily and is significantly similar to relaxin and relaxin-like factor. This gene is preferentially expressed in testis. Its expression in testis is restricted to interstitial cells surrounding seminiferous tubules, which suggests a role in sperm development and fertilization. [provided by RefSeq, Jul 2008],</p>