

IN35 Polyclonal Antibody

Cat No: HR1AP11823

For research use only

Overview

Product Name	IN35 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	31kDa
GeneID?Human?	IFI35 IFP35
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Not yet known.,induction:By interferon gamma.,similarity:Belongs to the NMI family.,subcellular location:Nuclear following IFN treatment.,subunit:Homodimer. Also interacts with B-ATF.,tissue specificity:In a wide range of cell types, including fibroblasts, macrophages, and epithelial cells.,