

IFNA2 Polyclonal Antibody

Cat No: HR1AP11816

For research use only

Overview

Product Name	IFNA2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	IFNA2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	interferon alpha 2(IFNA2) Homo sapiens This gene is a member of the alpha interferon gene cluster on chromosome 9. The encoded protein is a cytokine produced in response to viral infection. Use of the recombinant form of this protein has been shown to be effective in reducing the symptoms and duration of the common cold. [provided by RefSeq, Jun 2011],