

I20L2 Polyclonal Antibody

Cat No: HR1AP11813

For research use only

Overview

Product Name	I20L2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	ISG20L2 HSD-38 HSD38
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	interferon stimulated exonuclease gene 20 like 2(ISG20L2) Homo sapiens This gene encodes a 3'-5' exoribonuclease that may be involved in the processing of the 12S pre-rRNA. Pseudogenes have been identified on chromosomes 6 and 11. [provided by RefSeq, Dec 2014],