

# HSPB9 Polyclonal Antibody

Cat No: HR1AP11810

For research use only

## Overview

Product Name	HSPB9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GeneID?Human?	HSPB9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the small heat shock protein (HSP20) family.,tissue specificity:Testis specific.,