

## HSPB7 Polyclonal Antibody

Cat No: HR1AP11809

For research use only

### Overview

Product Name	HSPB7 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	18kDa
GeneID?Human?	HSPB7 CVHSP
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>similarity:Belongs to the small heat shock protein (HSP20) family.,subunit:Interacts with C-terminal domain of actin-binding protein 280.,tissue specificity:Isoform 1 is highly expressed in adult and fetal heart, skeletal muscle, and at a much lower levels in adipose tissue and in aorta. Undetectable in other tissues. Isoform 2 and isoform 3 are poorly detected in heart.,</p>