

## HS74L Polyclonal Antibody

Cat No: HR1AP11807

For research use only

### Overview

Product Name	HS74L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	HSPA4L APG1 OSP94
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>heat shock protein family A (Hsp70) member 4 like(HSPA4L) Homo sapiens The protein encoded by this gene is heat shock inducible and may act as a chaperone. The encoded protein can protect the heat-shocked cell against the harmful effects of aggregated proteins. This gene is highly expressed in leukemia cells and may be a good target for therapeutic intervention. Several transcripts encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2015],</p>