

HS12B Polyclonal Antibody

Cat No: HR1AP11806

For research use only

Overview

Product Name	HS12B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	75kDa
GeneID?Human?	HSPA12B C20orf60
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>heat shock protein family A (Hsp70) member 12B(HSPA12B) Homo sapiens The protein encoded by this gene contains an atypical heat shock protein 70 (Hsp70) ATPase domain and is therefore a distant member of the mammalian Hsp70 family. This gene may be involved in susceptibility to atherosclerosis. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2015],</p>