



HSPB3 Polyclonal Antibody

Cat No: HR1AP11808

For research use only

Overview

Product Name	HSPB3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	HSPB3 HSP27 HSPL27
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	heat shock protein family B (small) member 3(HSPB3) Homo sapiens This gene encodes a muscle-specific small heat shock protein. A mutation in this gene is the cause of autosomal dominant distal hereditary motor neuropathy type 2C.[provided by RefSeq, Sep 2010],