

## GRIP2 Polyclonal Antibody

Cat No: HR1AP11804

For research use only

### Overview

Product Name	GRIP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	114kDa
GeneID?Human?	GRIP2 KIAA1719
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>domain:PDZ 5 mediates the C-terminal binding of GRIA2 and GRIA3. PDZ 6 mediates interaction with the PDZ recognition motif of EFNB1. PDZ 7 mediates interaction with CSPG4.,function:May play a role as a localized scaffold for the assembly of a multiprotein signaling complex and as mediator of the trafficking of its binding partners at specific subcellular location in neurons.,similarity:Belongs to the GRIP2 family.,similarity:Contains 7 PDZ (DHR) domains.,subunit:Interacts with EFNB1, EFNB3, GRIA2, GRIA3, CSPG4 and GRIPAP1. Can form homomultimers and heteromultimers with GRIP1 (By similarity). Interacts with the C-terminal tail of PRLHR.,</p>