

HS12A Polyclonal Antibody

Cat No: HR1AP11805

For research use only

Overview

Product Name	HS12A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	74kDa
GeneID?Human?	HSPA12A KIAA0417
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the heat shock protein 70 family.,tissue specificity:Widely expressed with highest levels in brain, kidney and muscle.,