

GBP2 Polyclonal Antibody

Cat No: HR1AP11801

For research use only

Overview

Product Name	GBP2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	GBP2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>guanylate binding protein 2(GBP2) Homo sapiens This gene belongs to the guanine-binding protein (GBP) family, which includes interferon-induced proteins that can bind to guanine nucleotides (GMP, GDP and GTP). The encoded protein is a GTPase which hydrolyzes GTP, predominantly to GDP. The protein may play a role as a marker of squamous cell carcinomas. [provided by RefSeq, Jul 2013],</p>