

# FGFP1 Polyclonal Antibody

Cat No: HR1AP11800

For research use only

## Overview

|                       |   |
|-----------------------|---|
| Product Name          | FGFP1 Polyclonal Antibody   |
| Source                | Rabbit  |
| Applications          | WB,ELISA  |
| Species Reactivity    | Human   |
| Recommended Dilutions |   |
| Immunogen             |   |
| Species               | Rabbit  |
| Storage               | -20°C/1 year  |
| Isotype               |   |
| Clonality             |   |
| Concentration         | 1 mg/ml   |
| Observed band         | 25kDa   |
| GeneID?Human?         | FGFBP1 FGFBP HBP17  |
| Human Swiss-Prot No.  |   |
| Cellular localization |   |
| Alternative Names     |   |
| Background            | <p>fibroblast growth factor binding protein 1(FGFBP1) Homo sapiens This gene encodes a secreted fibroblast growth factor carrier protein. The encoded protein plays a critical role in cell proliferation, differentiation and migration by binding to fibroblast growth factors and potentiating their biological effects on target cells. The encoded protein may also play a role in tumor growth as an angiogenic switch molecule, and expression of this gene has been associated with several types of cancer including pancreatic and colorectal adenocarcinoma. A pseudogene of this gene is also located on the short arm of chromosome 4. [provided by RefSeq, Nov 2011],</p> |