

# FAIM1 Polyclonal Antibody

Cat No: HR1AP11797

For research use only

## Overview

Product Name	FAIM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	19kDa
GeneID?Human?	FAIM FAIM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	Fas apoptotic inhibitory molecule(FAIM) Homo sapiens The protein encoded by this gene protects against death receptor-triggered apoptosis and regulates B-cell signaling and differentiation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011],