

FBF1 Polyclonal Antibody

Cat No: HR1AP11798

For research use only

Overview

Product Name	FBF1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	124kDa
GeneID?Human?	FBF1 ALB KIAA1863
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,subunit:May interact with FAS cytoplasmic domain.,