

EMID2 Polyclonal Antibody

Cat No: HR1AP11796

For research use only

Overview

Product Name	EMID2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	EMID2 COL26A1 EMU2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	collagen type XXVI alpha 1 chain(COL26A1) Homo sapiens This gene encodes a protein containing an emilin domain and two collagen stretches. This gene may be associated with aspirin-intolerant asthma. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2013],