

DOT1L Polyclonal Antibody

Cat No: HR1AP11795

For research use only

Overview

Product Name	DOT1L Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	191kDa
GeneID?Human?	DOT1L KIAA1814 KMT4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	DOT1 like histone lysine methyltransferase(DOT1L) Homo sapiens The protein encoded by this gene is a histone methyltransferase that methylates lysine-79 of histone H3. It is inactive against free core histones, but shows significant histone methyltransferase activity against nucleosomes. [provided by RefSeq, Aug 2011],